

Proven Performance in Chiropractic Imaging

The Choice for Simplicity, Reliability and Value



AC1 17\" Chiropractic System Shown with Optional Electric Locks

**CHIROPRACTIC
RADIOGRAPHIC SYSTEMS**

AMERICOMP

Affordable Choices in Radiography

Since 1987, the quality standards of Americomp have been unsurpassed. Americomp systems remain the affordable choice in today's marketplace, offering the same long-term dependable performance that has become synonymous with the Americomp name. Our five-year manufacturer's parts warranty secures your investment today and guarantees your investment value tomorrow.

AC1 and AC2 Chiropractic System Features

Ideal for an existing practice or start-up practitioner, these reliable and versatile systems are an outstanding value.

The AC1 system accommodates up to a 14" x 17" (35 x 43 cm) cassette in either direction, allowing a broad range of procedures.

The AC2 system will accept up to a 14" x 36" cassette for full spine imaging. Both systems feature a traditional floor-wall or ceiling tubestand, and provide film focal distances of 40" to 72". The tubestand is constructed of heavy-duty steel, insuring many years of reliable service. Stainless steel ball bearings allow effortless positioning of the x-ray tube and precise alignment to the wallstand over the entire range of travel. The tubestand angulation dial is an easy-to-read, accurate indicator of tube angle. Sturdy locks provide secure placement of the x-ray tube. Optional electric locks allow quick tube positioning.

Both wallstands utilize a precision bearing system that provides effortless vertical travel of the grid and cassette tray. The full travel wallstands allow a wide range of exams including standing knee and ankle procedures. The wallstands feature a positive manual lock or optional electric lock.



AC2 36" Chiropractic System Shown with Optional Electric Locks



A3199 Mobile Radiographic / Adjusting Table Features

This versatile option to either x-ray system can be used as a chiropractic treatment table or as a mobile radiographic table. When used for adjustments or therapy, an optional 4" thick cushion with a facial slot provides superior patient comfort. Front wheel locks quickly secure the table into a stationary position. By removing the cushion, the table is ready for x-ray. It can be moved easily from room to room on its hospital-grade wheels that do not mar floors. The flat table top design has no obtrusive joints or fasteners and is available in two sizes; 24" x 56" or 24" x 76". Table heights of 20", 24" or 28" are available. A high quality, fine line grid and cassette tray are included.

High Performance Generators

Complement any Americomp radiographic system with a high-quality, state-of-the-art generator. Contact your Americomp dealer for details.

AMERICOMP

AMERICOMP INC.

Chicago, Illinois, USA

1-800-458-1782

Fax 773-583-1751

www.americompus.net

LITE 007-0505